



**RTX**




**PHILIPS**

**PHILIPS**

# RTX3466

Increase telemetry coverage with RTX3466 SMART HOPPING® Remote Antenna



Expanded coverage for  
challenging environments



Lower deployment costs



Simplified network  
infrastructure



Consistent performance  
in high-demand areas

## Reliable and secure wireless connectivity with WMTS in healthcare

Healthcare environments require robust and reliable wireless connectivity to telemetry devices, such as patient-worn devices, bedside monitors, infusion pumps, etc. This can be ensured over a reserved Federal Communications Commission (FCC) designated frequency band for US hospitals: Wireless Medical Telemetry Service (WMTS). WMTS is specifically allocated for transmitting patient medical data (e.g., vital signs) over wireless networks in healthcare environments.

This ensures that telemetry devices can reliably and securely transmit real-time data without interference from other wireless devices - vital for maintaining high-quality, uninterrupted patient care.

## Advancing WMTS with SMART HOPPING® for wireless patient monitoring

Although WMTS benefits from being a protected spectrum, incorporating the SMART HOPPING® technology further enhances its effectiveness by providing adaptive interference management, improved reliability, and low latency. SMART HOPPING® has specifically been developed to meet the unique demands of critical healthcare communication, making it an ideal solution for systems operating within the WMTS spectrum. Its reliability is especially crucial for telemetry devices, where uninterrupted communication is essential.

## Advantages of SMART HOPPING® in wireless patient monitoring systems:



### Interference avoidance

SMART HOPPING® adaptively switches frequencies and time slots upon detecting interference, ensuring a robust signal, and minimizing data loss.



### Efficient spectrum utilization

By channelizing signals into time slots, SMART HOPPING® allows multiple devices to share the same spectrum efficiently, allowing up to 1,024 devices on a single network.



### Deterministic roaming

SMART HOPPING® assesses signal strength, packet error rates, and fading conditions to determine optimal roaming times and locations, enhancing connectivity and sustaining low latency.



### Data confidentiality

Using a proprietary protocol, SMART HOPPING® avoids transmitting patient identifiers over the wireless interface.



### Bi-directional communication

SMART HOPPING® features bi-directional communication to also enable features such as “Find Device,” “Device Location,” and “Control Clinical Measurements” to enhance functionality of telemetry devices.

These advantages collectively contribute to a reliable, scalable, and secure wireless monitoring environment in healthcare settings.

## Extend wireless coverage with the SMART HOPPING® Remote Antenna

**RTX3466 – SMART HOPPING® Remote Antenna 1.4GHz** utilizes SMART HOPPING® technology to significantly boost the wireless range and connectivity of the **RTX3481 – SMART HOPPING® Access Point 1.4GHz**. Ideal for environments requiring broad and dependable wireless access, the RTX3481 Access Point can connect with up to two RTX3466 Remote Antennas to enlarge its coverage zone. Each RTX3466 Remote Antenna can be installed up to 68 feet (20 meters) away, extending radio signals for optimized performance. This setup delivers dependable, scalable wireless communication - critical for supporting healthcare staff, devices, and patients.

For detailed specifications on the RTX3481 Access Point, consult the product sheet.

### Expanded coverage for challenging environments

Hospitals are often difficult for wireless networks due to their structural complexity - multi-story buildings, reinforced concrete, shielded rooms like X-ray areas, and a variety of medical equipment that obstructs standard Wi-Fi signals. The RTX3466 Remote Antenna extends the reach of an access point into areas that would otherwise remain unserved, eliminating dead zones in essential spaces such as ICUs, operating theaters, distant wards, and restrooms.

By bridging physical barriers, the RTX3466 Remote Antenna ensures dependable connectivity for both clinical workflows and hospital operations.

### Lower deployment costs

In scenarios where the RTX3466 Remote Antenna is not used, achieving full coverage typically requires deploying additional access points - each demanding cabling and installation efforts. Strategically placed RTX3466 Remote Antenna(s) in the core switching infrastructure allows a single access point to cover a broader area, reducing the need for extra equipment.

This streamlined approach minimizes expenses related to hardware, wiring, and maintenance.

### Simplified network infrastructure

A dense deployment of access points to reach difficult areas can create signal overlap, interference, and complex network management. The RTX3466 Remote Antenna enables more targeted coverage, allowing access points to remain in convenient, accessible areas (such as server rooms or above ceilings) while the antennas are extended to where signal is most needed.

This results in a cleaner network design that is easier to troubleshoot and maintain - especially vital in healthcare settings where reliability is non-negotiable.

### Consistent performance in high-demand areas

Hospitals depend on wireless infrastructure to support a range of critical functions - from Electronic Health Record (EHR) systems and VoIP to telemetry and patient entertainment platforms. The RTX3466 Remote Antenna helps balance user load, minimizes network congestion, and ensures strong signal delivery, even in crowded or high-traffic zones.

It guarantees high-capacity, stable wireless performance for vital medical applications and simultaneous users.



# Technical specifications

Physical	
Chassis (only dimensions (H x W x L))	<ul style="list-style-type: none"> <li>30 mm x 163 mm x 163 mm (1.2 in x 6.38 in x 6.38 in)</li> </ul>
Weight with internal antennas	<ul style="list-style-type: none"> <li>&lt;~300 g (~0.66 lb)</li> </ul>
Mounting options	<ul style="list-style-type: none"> <li>Wall mount (standard)</li> <li>Wall mount with cosmetic cover</li> <li>Flush ceiling mount with cosmetic ring</li> <li>Mount below ceiling tile</li> <li>Mount below ceiling tile with cosmetic cover</li> </ul>
Environmental	
Operating temperature	<ul style="list-style-type: none"> <li>0 to +55°C (32 to 131°F)</li> </ul>
Storage temperature	<ul style="list-style-type: none"> <li>-40 to +70°C (-40 to 158°F)</li> </ul>
Humidity range (operating)	<ul style="list-style-type: none"> <li>15% RH to &lt; 95% RH @ 40°C non-condensing</li> </ul>
Humidity range (storage)	<ul style="list-style-type: none"> <li>5% RH to &lt; 90% RH @ 60°C non-condensing</li> </ul>
Altitude	<ul style="list-style-type: none"> <li>-500 meters (-1,640 ft) to 3,048 meters (10,000 ft) (operating)</li> <li>-500 meters (-1,640 ft) to 4,700 meters (-15,420 ft) (non-operating) (106 Kpa to 56 Kpa)</li> </ul>
Electrical	
Power (input)	<ul style="list-style-type: none"> <li>Power input 5V nom (4.5 - 5.5) 110mA (from RTX3481 Access Point)</li> </ul>
LED indicators	Two LEDs for: <ul style="list-style-type: none"> <li>Power</li> <li>Radio activity</li> </ul>
SMART HOPPING® radio	
SMART HOPPING® Remote Antenna frequency range	<ul style="list-style-type: none"> <li>1390 - 1395 MHz + 1395 - 1400 MHz + 1427 - 1432 MHz + 1432 - 1435 MHz</li> <li>WMTS - Wireless Medical Telemetry Service (USA)</li> </ul>
Antenna type	<ul style="list-style-type: none"> <li>SMART HOPPING® dual internal antennas</li> </ul>
Space diversity	<ul style="list-style-type: none"> <li>Uses two antennas and selects the antenna with the best signal</li> </ul>
Frequency diversity	<ul style="list-style-type: none"> <li>Dynamically selects RF channels for best signals</li> </ul>
Compatibility mode	
Channel spacing	<ul style="list-style-type: none"> <li>1.6 MHz</li> </ul>
Modulation	<ul style="list-style-type: none"> <li>GFSK</li> </ul>
Power output	<ul style="list-style-type: none"> <li>12.5 (±3.0) dBm (12.6 mW to 25 mW)</li> </ul>
Antenna gain	<ul style="list-style-type: none"> <li>2 dBi</li> </ul>
Time slot types supported	<ul style="list-style-type: none"> <li>Single or long</li> </ul>
Total time slots available	<ul style="list-style-type: none"> <li>One RTX3466 connected: 9 (and 9 for the Remote Antenna)</li> <li>Two RTX3466 connected: 6 (and 6 for each Remote Antenna)</li> <li>Time slots allocated for roaming: 14</li> </ul>
RF sensitivity @ 10-3 BER	<ul style="list-style-type: none"> <li>Min. -87 dBm (nom. -92 for RTX3481 Access Point)</li> </ul>
RSSI level per cell radius	<ul style="list-style-type: none"> <li>-67 dBm</li> </ul>
Advanced mode (4 channel mode)	
Channel spacing	<ul style="list-style-type: none"> <li>1.728 MHz</li> </ul>
Modulation	<ul style="list-style-type: none"> <li>PI/2-DBPSK, PI/4-DQPSK</li> </ul>
Power output	<ul style="list-style-type: none"> <li>High: 14 (±3) dBm in RTX3481 Access Point</li> </ul>
Antenna gain	<ul style="list-style-type: none"> <li>2 dBi</li> </ul>

# Technical specifications

## Advanced mode (4 channel mode) (continued)

Time slot types supported	<ul style="list-style-type: none"> <li>• Single or long</li> </ul>
Time slots allocated for wireless clients	<ul style="list-style-type: none"> <li>• One RTX3466 connected: 9 (and 9 for the Remote Antenna)</li> <li>• Two RTX3466 connected: 6 (and 6 for each Remote Antenna)</li> <li>• Time slots allocated for roaming: 14</li> </ul>
RF sensitivity @ 10-3 BER	<ul style="list-style-type: none"> <li>• Min. -87 dBm for BPSK</li> <li>• Min. -86 dBm for QPSK</li> <li>• Nom. -92 for RTX3481 Access Point</li> </ul>

## Electrical installation

Fire safety	<ul style="list-style-type: none"> <li>• RTX3466 – SMART HOPPING® Remote Antenna 1.4 GHz is listed for use within “Other Spaces Used for Environmental Air (Plenum)” per NFPA70: 2011, Article 300.22. Note: The term ‘plenum’ as used in Article 300.22 Section C correlates with the use of the term ‘plenum’ in NFPA 90A-2009, Standard for the Installation of Air-Conditioning and Ventilating Systems, and other mechanical codes where the plenum is used for return air purposes, as well as some other air-handling spaces. The area above dropped ceilings is an example of plenum space.</li> </ul>
-------------	--

## Interface connections

Control signals	<ul style="list-style-type: none"> <li>• RJ-45 female socket</li> </ul>
Core AP 75 Ohm coaxial cable connector	<ul style="list-style-type: none"> <li>• One standard 75 Ohm connector is provided for the coaxial cable on which RF and DC sense signals are carried to the RTX3466 Remote Antenna.</li> </ul>

## Ordering details

RTX3466 - SMART HOPPING® Remote Antenna 1.4 GHz	<ul style="list-style-type: none"> <li>• Contact RTX for more details at <a href="mailto:sales@rtx.dk">sales@rtx.dk</a></li> </ul>
Part number	<ul style="list-style-type: none"> <li>• RTX3466 SH2.0 Remote Antenna 1.4 GHz: 94100010</li> <li>• Philips Smart-hopping 2.0 Remote Antenna 1.4 GHz OEM: 989803221641</li> </ul>

# Application Standards

## Radio as accessory to RTX3481 (ITS867216A) Standard

FCC Part 95H

## As accessory to ITS4843E Standard

FCC Part 95H

## EMC as accessory to RTX3481 (ITS867216A) Standard

EN 55032:2015+A11:2020+A1:2020

EN 55035:2017+A11:2020

EN IEC 61000-3-2:2019+A1:2021

EN 61000-3-3:2013+A1:2019+A2:2021

BS EN 60601-1-2:2015+A1:2021

FCC 15 Subpart B

## EMC as accessory to ITS4843E Standard

EN 55032:2015+A11:2020+A1:2020

EN 55035:2017+A11:2020

EN IEC 61000-3-2:2019+A1:2021

EN 61000-3-3:2013+A1:2019+A2:2021

BS EN 60601-1-2:2015+A1:2021

FCC 15 Subpart B

## Safety Standard

### Area

CAN/CSA C22.2 No. 62368-1:19  
UL 62368-1:Ed.3,2019(R2021)

- Safety

FCC 47 CFR §2.1091

- SAR

UL 2043, FIFTH EDITION

- Smoke Release